

MezoPharma PEEL

About MezoPharma acids:

The innovation of MezoPharma PEEL acids is the use of a specialized base based on plant glycol, which allows even more effective MezoPharma peels.

MezoPharma Glycolic Acid 30% pH 1.0	
Description:	Glycolic acid is recommended for skin of all ages. Especially for dry skins and requiring regeneration, because it moisturizes by retaining water in the epidermis. It improves the elasticity and firmness of the skin, while stimulating it to rebuild.
For whom:	For people with gray skin that requires moisturizing and regeneration; for people with a photoaging problem and excessive wrinkles.
When to use:	Recommended as preparation for rejuvenating therapies and for combination with aesthetic treatments
Capacity:	50ml pipette bottle
Application time:	3-5 minutes

MezoPharma Glycolic Acid 50% pH 0.8	
Description:	Glycolic acid is recommended for skin of all ages. Especially for dry skins and requiring regeneration, because it moisturizes by retaining water in the epidermis. It improves the elasticity and firmness of the skin, while stimulating it to rebuild.
For whom:	For people with a photoaging problem and excessive wrinkles.
When to use:	Recommended as intensive rejuvenation therapies
Capacity:	50ml pipette bottle
Application time:	2-4 minutes

MezoPharma Lactic Acid 30% pH 1.6	
Description:	Lactic acid belongs to the group of alpha-hydroxy acids and is obtained by fermenting sugar by lactic acid bacteria (<i>Lactobacillus</i>).
For whom:	For people with hyperpigmentation, gray and earthy skin, smoker's skin, for combined therapies
When to use:	Recommended as depigmentation and brightening therapies
Capacity:	50ml pipette bottle
Application time:	3-8 minutes

MezoPharma Lactic Acid 50% pH 1.4	
Description:	Lactic acid belongs to the group of alpha-hydroxy acids and is obtained by fermenting sugar by lactic acid bacteria (<i>Lactobacillus</i>).
For whom:	For people with hyperpigmentation, gray and earthy skin, smoker's skin,
When to use:	Recommended as depigmentation and brightening therapies
Capacity:	50ml pipette bottle
Application time:	3-6 minutes

MezoPharma Azeloglycine 10% pH 5.0

Description:	Azeloglycine is a derivative of azelaic acid in the form of a liquid solution, easily water soluble. It was created by combining azelaic acid with an amino acid glycine. Such a modification allows obtaining an active substance with better physical and chemical properties.
For whom:	For people with discoloration in any form, for skin requiring regeneration and color harmonization
When to use:	Recommended as brightening, refreshing therapies, as well as moisturizing and protecting particularly sensitive skin.
Capacity:	50ml pipette bottle
Application time:	10-20 minutes
<i>Safe for pregnant women</i>	

MezoPharma Pyruvic Acid 40% pH 1.2

Description:	Pyruvic acid is used for skin that requires intensive regeneration and repair. It works bacteriostatic, regulates sebum secretion and reduces adherence between cells in the epidermal layer of the epidermis, which in turn causes intensive peeling of the epidermis.
For whom:	For people with all forms of acne, for people with hyperkeratosis, scars, acne, extensive skin discoloration.
Capacity:	50ml pipette bottle
Application time:	1-4 minutes

MezoPharma Ferulic Complex pH 1.8

Description:	Ferulic Complex is an innovative combination of 5% ferulic acid, 10% lactic acid and 6% glycolic acid.
For whom:	For mature skin that needs refreshing and lifting.
When to use:	Complex dedicated to intensive rejuvenation and lifting treatments.
Capacity:	50ml pipette bottle
Application time:	3-6 minutes

MezoPharma Pyruvic Complex pH 1.8

Description:	Pyruvic Complex is a combination of 20% pyruvic acid and 10% lactic acid.
For whom:	Complex dedicated to impure, oily, uneven skin and excessive work of sebaceous glands or as intensive dermoplasty of demanding skin rejuvenation.
When to use:	Should be used as intensive depigmentation or anti-acne therapies.
Capacity:	50ml pipette bottle
Application time:	4-8 minutes

MezoPharma Neutralizer

Description:	The high content of aloe and d-panthenol soothes and rebuilds skin functions.
When to use:	Neutralizing and soothing preparation dedicated as a pH restoration after treatments acid.
Capacity:	250ml bottle
Application time:	10 minutes

Mezopharma *Hyal* is an innovative line of advanced mesotherapy cocktails. Each and every one of its products is an advanced meso cocktail rich in biomimetic peptides, hyaluronic acid, aminoacids and vitamins.


Each cocktail is individually sterilised in an advanced processes equal to those used in pharmacy. Sterilising MezoPharma Hyal cocktails is a joint process of aseptic filtration and heat sterilization. It is one of the most innovative sterilising processes in biotechnological industry which maintains full effectiveness of the heat-sensitive substances while maintaining maximal safety of your client.

Safety



Modern associated strilisation system

Advanced Receptures



Biomimetic peptides, one of the newest biotechnological discoveries

Packaging



Single-use and sterile vials protected from detrimental UV effects.

	AQUALIFT COCKTAIL	
	Aqualift is a cocktail for face and neck mesotherapy that effectively prevents first signs of aging in the form of mimic wrinkles. Hyaluronic acid included in Aqualift cocktail fully recovers the water barrier moisturizing and brightening the skin.	
	Frequency	1 treatment per 14 days
	Application area	face, neck, décolletage
	Results	moisturizing, lifting, improvement of the skin's condition
	Target group	first signs of ageing , dry skin, loose skin, recommended for age group 30+
Precautions	combination with other cocktails in one treatment is forbidden	

	BOTOLIFT COCKTAIL	
	Botolift cocktail is an exclusive mesotherapy preparation dedicated to skin rejuvenation and firming . Numerous peptides and biomimetic peptides relax the muscles responsible for deep creases and wrinkles. Botolift cocktail contains advanced forms of amino acids and coenzymes which delay all aging processes.	
	Frequency	1 treatment per 14 days
	Application area	face, neck, décolletage
	Results	relaxes muscles responsible for deep lines and wrinkles; Acts as botulinum toxin but does not demonstrate high risk of adverse effects (muscle damage, prolapse, or production of antibodies)
	Target group	addressed for age group 50+ for advanced skin rejuvenation and lifting
Precautions	combination with other cocktails in one treatment is forbidden	



HAIR COCKTAIL

Advanced mesotherapy cocktail for the head skin used to **eliminate baldness and excessive hair loss**. The innovation of Hair Cocktail consists of including 6 biomimetic peptides, more than 28 aminoacids, many vitamins and innovative troxerutin and trifolium pratense flower extract.

Frequency	1 treatment per 7-9 days
Application area	head and eyebrow skin
Results	reduced hair loss, reduced baldness, elimination of seborrheic dermatitis on the scalp
Target group	people suffering from baldness, excessive hair loss, increased sebaceous gland activity over the scalp
Precautions	none



STRETCH COCKTAIL

Stretch cocktail is dedicated to neck and body mesotherapy. The dimethylaminoethanol tartare contained in Stretch Cocktail prevents tissue and skin damage. **It reduces the appearance of stretch marks and tightens the skin**. It positively affects the fibroblast metabolism inducing collagen and elastin synthesis.

Frequency	1 treatment per 14-21 days
Application area	wide stretch marks all over the body, loose skin, scars
Results	elimination of stretch marks, scar softening, improved skin density and elasticity
Target group	people with loose skin, stretch marks and scars
Precautions	none



VITALIFT COCKTAIL

Vitalift is a cocktail for face, neck and neckline mesotherapy. **It improves and evens out the skin tone** from the first application. The advanced formula of the cocktail removes discolorations and inhibits the melanogenesis process preventing new discolorations.

Frequency	1 treatment per 14-21 days
Application area	loose skin, stretch marks and scars
Results	elimination of stretch marks, scar softening, improved skin density and elasticity
Target group	dedicated for people with skin's discoloration, uneven skin tone, acne skin, smoker's and oily skins
Precautions	combination with other cocktails in one treatment is forbidden.



With the massive advance in technology and biotech industry we currently observe a substantial advantage of using High Tech in aesthetic procedures. This is caused by very high effectiveness of the procedures as well as minimal risk of any kind of complications.

To meet our customers needs a specialistic OXY line was created based on a hyaluronic acid with exceptionally low molecular weight – LMHA (low mollecular hyaluronioc acid). Products that create MezoPharma OXY line are innovative and sterile cocktails dedicated especially to High Tech procedures such as oxygen infusion, radio-frequency, radio-frequency micro needle resurfacing, sonophoresis, electroporation etc.


Commonly used hyaluronic acid does not have and transdermal qualities to itself – which is caused by too large molecules used in the production process. MezoPharma OXY line of products contain the low molecular hyaluronic acid that had been specifically selected to perform better. It has substantially better ability to penetrate deeper layers of skin in comparison to standard components. This way LMHA stimulates collagen and elastin synthesis, multiplies fibroblats and envigorates all rejuvenating processes in a way faster and more effective manner.

Safety



Modern associated strilisation system

Advanced Receptures



Biomimetic peptides, one of the newest biotechnological discoveries

Packaging



Single-use and sterile vials protected from detrimental UV effects.



HydraVIT

Perfectly eliminates all the first processes and signs of skin aging; due to the high content of vitamin C it strengthens the capillaries and directly brightens discoloration.

Frequency	1 treatment per 12-14 days
Composition	Vitamin C 10%, 5% hyaluronic acid, 5% retinol, vitamin E, mannitol, dextran, sodium chloride
Target group	for clients with the first signs of aging, for smokers' skins, perfect for eye revitalization, for acne skin
Precautions	do not combine with other cocktails in one treatment



HydraLIFT

Cocktail with a strong firming effect, improves skin density and elasticity, eliminates surface wrinkles associated with facial expressions

Frequency	1 treatment per 12-14 days
Composition	Hyaluronic acid 2.5% l.m., methylsilanol mannuronate, glutathione, DMAE 2%, sodium chloride, mannitol
Target group	for people who require strong skin firming, for people over 40 years of age, for dry and mature skin
Precautions	do not combine with other cocktails in one treatment



HydraDERM

Relaxes the muscles responsible for deep furrows and wrinkles, procedure compatible with botulinum toxin, but without the risk of complications (muscle damage, falling or producing antibodies)

Frequency	1 treatment per 12-14 days
Composition	3.5%, hyaluronic acid, 9 biomimetic peptides, 16 amino acids, complex of B vitamins, DMAE 2%, simple peptides, numerous coenzymes
Target group	an exclusive formulation for very demanding customers, dedicated to people over 50 as a deep regeneration and biomimetic skin lifting
Precautions	do not combine with other cocktails in one treatment